

Nuvo⁸ Mark 5



Nuvo¹⁰ Mark 5

Nuvo 8

115 Volt **230 Volt**
60 Hz **50 Hz**

580 **980**

585 **985**

Standard Model Part Number

OCSI Model Part Number
(Oxygen Concentration Status Indicator)

Nuvo 10

230 Volt
50 Hz

1000

1005

625 Watts 490 Watts

90% (+5% / -3%)
at 8 LPM

2 to 8 LPM

15 PSIG (1.2 BAR)

55.5 lbs / 25.2 kg

16"W x 16"D x 28"H
39 cm x 40 cm x 71 cm

Average Power Consumption

Oxygen Concentration

Liter Flow

Outlet Pressure

Weight

Dimensions

600 Watts

90% (+5% / -3%)
at 10 LPM

2 to 10 LPM

20 PSIG (1.38 BAR)

64.5 lbs / 29.3 kg

16"W x 16"D x 30"H
39 cm x 40 cm x 76 cm

Accessories available upon request -- Please contact your sales rep!

Additional Features



- Built on the Reliable Nuvo Platform
- Easy Access to filter and hour meter
- Superior Grade of Molecular Sieve
- Spring Mounted High Efficiency Compressor

2010-8406 REV A