

# Nuvo<sup>Mark 5</sup> *Lite<sup>3</sup>*



# Nuvo<sup>Mark 5</sup> *Lite*

## Nuvo Lite 3

230 Volt

50 Hz

--

935

Standard Model Part Number

OCSI Model Part Number  
(Oxygen Concentration Status Indicator)

## Nuvo Lite

115 Volt

60 Hz

230 Volt

50 Hz

520

525

920

925

180 Watts

90% (+6.5% / -3%)  
at 3 LPM

0.125 to 3 LPM

40 dBA typical

7.0 PSI

32 lbs / 14.5 kg

14"W x 9"D x 23"H  
36 cm x 23 cm x 58.5 cm

### Average Power Consumption

### Oxygen Concentration

Liter Flow

Sound Level

Outlet Pressure

Weight

Dimensions

330 Watts

90% (+6.5% / -3%)  
at 5 LPM

0.125 to 5 LPM

45 dBA typical

7.0 PSI

32 lbs / 14.5 kg

14"W x 9"D x 23"H  
36 cm x 23 cm x 58.5 cm

300 Watts

**Accessories available upon request -- Please contact your sales rep!**

## Additional Features



Lockable Flow Control Valve



Metal Outlet Fitting



Resettable Circuit Breaker

**NIDEK**  
*Medical*

CE  
0413