

EASYPULSE *TOC*

Total Oxygen Concentrator



The Perfect Balance of Portable and Stationary Oxygen Delivery

Total Oxygen Concentrator

It is the perfect balance of mobile and stationary oxygen delivery. It delivers continuous flow and pulse flow all in the same unit weighing less than 11.5 pounds. No need to have 2 separate oxygen concentrators.

Easy to Use

The EasyPulseTOC has an easy to use patient interface. One simple knob easily sets the oxygen flow to either pulse or continuous mode. Unlike many other TOC's, there are no confusing electronic controls.

Perfect for the Non-Delivery Oxygen Model

The EasyPulseTOC is the ideal product to help you reduce costly cylinder deliveries and eliminate the need for stationary oxygen concentrators.

Makes Mobility Easy

The convenient fold-up cart with its adjustable telescoping handle and oversized wheels is easily attached to the TOC enabling the patient to take the unit anywhere and by any means. It is FAA approved for use on airlines and fits easily under an airline seat.

All You Need is in the Box

Open the box and you are set to go! It contains the Concentrator, Cart, AC Power Supply, AC Power Cord, DC Power Cord, Battery, Accessory Bag and User Manual.

Excellent Battery Life

The battery lasts 3 hours. And with optional, easy to change extra batteries, you can extend your mobility as long as you wish.

Made in



USA



Easy to Use Patient Interface

SPECIFICATIONS

Dimensions:

Height 14.63 in. (37.16 cm)
 Width 10.50 in. (26.67 cm)
 Depth 7.00 in. (17.78 cm)

Weight:

11.5 lbs (5.2 kg)

Operating Temperature:

TOC: 41° to 104°F (5° to 40°C)
 Charging Battery: 32° to 113°F (0° to 45°C)

Altitude:

Sea Level to 10,000 ft (3048 m)

Specifications are subject to change without prior notice.

ORDERING INFORMATION

Part Number	Description
PM4400	EasyPulse TOC Oxygen Concentrator (includes Concentrator, Cart, AC Power Supply, AC Power Cord, DC Power Cord, Battery, Accessory Bag and User Manual)
508291	AC Power Cord
508188	AC Power Supply
508189	DC Power Cord
508215	Battery
508260	Cart
508301	Accessory Bag
508306	Desktop Charger Kit



**ISO 13485
 Certified**

Precision Medical, Inc.
 300 Held Drive • Northampton, PA 18067 USA
 Phone: (+001) 610-262-6090 Toll-Free Phone: 1-800-272-7285
 Fax: (+001) 610-262-6080 Toll-Free Fax: 1-800-352-1240
www.precisionmedical.com